



 **ORFACE**
ORTHOPEDIC PRODUCTS **NECK COLLAR GROUP**



SPONGE NECK COLLAR

SÜNGER BOYUNLUK
ВОРОТНИК С ГУБЧАТОЙ ШЕЕЙ

OB 300



Provides superior comfort and support to the neck with its anatomically designed structure. Limits excessive cervical movement and provides additional warmth around the neck. MR compatible.

INDICATIONS

- Cervical disc hernia
- Muscle injuries
- Cervical arthrosis
- Muscle spasms
- Post-traumatic degenerative changes of the cervical spine



XS	S	M	L
22-27	28-33	34-39	40-45



SPONGE NECK COLLAR WITH UNDERWIRE SUPPORT

(SÜNGER BOYUNLUK BALEN DESTEKLİ)
ВОРОТНИК С ГУБЧАТОЙ ШЕЕЙ И ПОДДЕРЖКОЙ НА КОСТОЧКАХ

OB 301



It provides superior comfort and support to the neck through anatomical designed. It refrains from excessive cervical activity and provides extra warmth around the neck. Compatibility MR.

INDICATIONS

- Cervical disc hernia
- Muscle injuries
- Cervical arthrosis
- Muscle spasms
- Post-traumatic degenerative changes of the cervical spine

XS	S	M	L
22-27	28-33	34-39	40-45





NELSON NECK COLLAR

NELSON BOYUNLUK
ШНУРОК НЕЛЬСОНА

OB 303



It provides superior comfort and support to the neck through anatomical designed. It helps the patient during the treatment by restricting the movements in the neck area. Compatibility MR.

INDICATIONS

- Cervical disc hernia
- Muscle injuries
- Cervical arthrosis
- Muscle spasms
- Post-traumatic degenerative changes of the cervical spine

XS	S	M	L
22-27	28-33	34-39	40-45



NELSON NECK COLLAR COVERED FABRIC

NELSON BOYUNLUK KUMAŞ KAPLI
БОРОТНИК НЕЛЬСОНА, ОБЯНУТЫЙ ТКАНЬЮ

OB 304



It provides superior comfort and support to the neck through anatomical designed. It helps the patient during the treatment by restricting the movements in the neck area. Compatibility MR.

INDICATIONS

- Cervical disc hernia
- Muscle injuries
- Cervical arthrosis
- Muscle spasms
- Post-traumatic degenerative changes of the cervical spine

S	M	L	XL
32-35	36-39	40-43	43-46





PHILADELPHIA NECK COLLAR

PHILADELPHIA BOYUNLUK
 ШЕЙНЫЙ БАНДАЖ «ФИЛАДЕЛЬФИЯ»

OB 306



It is a 2-piece neck collar. It helps the patient during the treatment by limiting flexion, lateral and rotational movements of the neck. It has a chin support. Compatibility MR.

INDICATIONS

- Soft tissue traumas
- Stable bone and ligament injuries
- Post-trauma, recovery and preservation purposes
- Cervical disc irritation
- Suspected cervical trauma in unconscious patients

S 32-37
M 38-42
L 42-46



PHILADELPHIA NECK COLLAR TRACHEOSTOMY OPEN

PHILADELPHIA BOYUNLUK
 TRAKEOSTOMISI AÇIK

ОТКРЫТАЯ ТРАХЕОСТОМИЯ ШЕИ В ФИЛАДЕЛЬФИИ

OB 307



It is a 2-piece neck collar. It helps the patient during the treatment by limiting flexion, lateral and rotational movements of the neck. It is a collar with trachea open. Compatibility MR.

INDICATIONS

- Soft tissue traumas
- Stable bone and ligament injuries
- Post-trauma, recovery and preservation purposes
- Cervical disc irritation
- Suspected cervical trauma in unconscious patients

S 32-37
M 38-42
L 42-46





CERVICAL CHIN NECK COLLAR

SERVİKAL ÇENELİ BOYUNLUK
 ШЕЙНЫЙ БАНДАЖ ШЕЙНОЙ ЧЕЛЮСТИ

OB 308



Anatomically designed to the neck form and manufactured from plastic material. Height is adjustable with velcro fasteners. Stabilizes the neck area and relieves the pressure on the neck muscles.

INDICATIONS

- Cervical disc hernia
- Cervical arthrosis
- Post-trauma, recovery and preservation purposes

S	M	L	XL	2XL
32-37	36-39	40-43	43-46	46-50



WIRE FRAME NECK COLLAR

WIRE FRAME BOYUNLUK
 КАРКАСНЫЙ ШНУРОК ИЗ КЕРАМИКИ

OB 310



It is a 2-piece neck collar. It helps the patient during the treatment by limiting flexion, lateral and rotational movements of the neck. It is a collar with trachea open. The back of the neck is supported by a sponge fabric.

INDICATIONS

- Cervical disc hernia
- Cervical arthrosis
- Post-trauma, recovery and preservation purposes

S	M
33-39	40-46





EMERGENCY NECK COLLAR

ACİL YARDIM BOYUNLUĞU
АВАРИЙНЫЙ ШНУРОК

OB 311



It is designed to use for emergency patient. Manufactured from plastic material. Inner side is coated by soft foam material. 3 levels of height adjustment possible with security lock system and trachea open. Compatibility MR.

INDICATIONS

- For first aid and stabilization in emergencies

S	M
20-46	28-58



MINEVRA ORTHOSIS

MİNERVA ORTEZİ
МИНЕРВА ОРТЕЗИ

OB 312



It is a 2-pieces neck collar and body orthosis. The inner side is covered with a sponge terry fabric and the outer surface is made of polyethnylene material and is supported with aluminum. It fixed the chest area and restricts the extension, flexion, lateral an rotation movements of the neck.

INDICATIONS

- Cervical spine trauma
- Vertebral fractures
- Cervical arthrosis
- Post-operative rehabilitation
- Degenerative joint diseases

S	M	L
33-37	38-42	42-46
70-90	90-110	110-130





MINEVRA OSTOSIS SHORT

KISA MİNERVA ORTEZİ
 КИСА МИНЕРВА ОРТЕЗ

OB 314



It is 2-pieces neck collar and body orthosis. The inside is covered with a sponge terry fabric and the outer surface is made of polyethylene material and is supported with aluminum. It fixes the chest area and restricts the extension, flexion, lateral flexion and rotation movements of the neck.

INDICATIONS

- Cervical spine trauma
- Vertebral fractures
- Cervical arthrosis
- Post-operative rehabilitation
- Degenerative joint diseases

S	M	L
32-37	38-42	42-46



MIAMI NECK COLLAR

MİAMİ BOYUNLUK
 ШЕЙНЫЙ ВОРОТНИК МАЙАМИ

OST 1001



It is made of 2 plastic parts and has sponge supports covered with fabric on the inside. Thanks to the 6 different height levels on the top, height adjustment can be made between 9 cm and 16 cm. The tracheostomy hole is an open neck collar.

INDICATIONS

- Sprain
- Postoperative immobilization
- Soft tissue injuries
- Severe cervical syndrome
- Relaxation of the vertebral structure

UNIVERSAL SIZE

