



ORFACE

ORTHOPEDIC ^{FT} SHOULDER ARM SUPPORT



EPICONDYLITIS PRACTICAL SYSTEM

ΕΠΙΚΟΝΔΙΛΙΤ ΒΑΝΔΙ ΠΡΑΤΙΚ ΣΙΣΤΕΜΛΙ
 ЭПИКОНДИЛИТ ШИНА ПРАКТИЧЕСКАЯ СИСТЕМА

OD 200



It is made of leather and breathable fabric. With special 3-points pad that provides local coompression, reduce ache on the tendon. It provides easy usage with practical lacing mechanism. It is designed as suitable for use on both sides. It accelerates blood circulation and maintains body temperature.

INDICATIONS

- Tennis elbow
- Golfer's elbow

S
22-26

M
26-32



EPICONDYLITIS BANDAGE

ΕΠΙΚΟΝΔΙΛΙΤ ΒΑΝΔΑΛΙ
 БАНДАЖ ОТ ЭПИКОНДИЛИТА

OST 1002



It provides local compression and reduce ache on the tendon with special pad. It is designed as suitable for use on both sides. It accelerates blood circulation and maintains body temperature.

INDICATIONS

- Golfer's elbow
- Tennis elbow

UNIVERSAL SIZE





EPICONDYLITIS SUPPORTED ELBOW SUPPORT

EPİKONDİLİT DESTEKLİ DİRSEKLİK
 ПОДДЕРЖКА ЛОКОТЬЯ С ПОДДЕРЖКОЙ
 ЭПИКОНДИЛИТА

OST 1003



With special 3-points pad that provides local coompression, reduce ache on the tendon. It is designed as suitable for use on both sides. It accelerates blood circulation and maintains body temperature.

INDICATIONS

- Golfer's elbow
- Tennis elbow



UNIVERSAL SIZE



EPICONDYLITIS SUPPORTED ELBOW SUPPORT

EPİKONDİLİT DESTEKLİ DİRSEKLİK
 ПОДДЕРЖКА ЛОКОТЬЯ С
 ПОДДЕРЖКОЙ ЭПИКОНДИЛИТА

OST 1004



With special 3-points pad that provides local coompression, reduce ache on the tendon. It is designed as suitable for use on both sides. It accelerates blood circulation and maintains body temperature.

INDICATIONS

- Golfer's elbow
- Tennis elbow

UNIVERSAL SIZE





TENNIS ELBOW

TENİSÇİ DİRSEKLİĞİ
ТЕННИС ЭЛББОУ

OD 201



It is made of neoprene, via its special silicone pad that provides local compression, it reduces the pain in the tendon. Special design allows for both left and right arm use, it accelerates blood circulation and preserves body temperature.

INDICATIONS

- Lateral Epicondylitis
- Golfer's Elbow
- Tennis player's elbow
- In painful elbow conditions

S	M	L	XL	2XL
22-24	24-26	26-28	28-30	30-32



KNITTED SILICONE SUPPORTED ELBOW

ÖRME SİLİKON DESTEKLİ DİRSEKLİK
ТРИКОТАЖНЫЙ ЛОКОТЬ
С СИЛИКОНОВОЙ ПОДЛОЖКОЙ

OD 202



Manufactured from respirable and elastic woven fabric and designed anotomically. It has droplet featured silicone pads that provides compression on both sides of elbow. Local compression reduces tendon tension and so the pain.

INDICATIONS

- Golfer's elbow
- Tennis elbow

S	M	L	XL	2XL
22-24	24-26	26-28	28-30	30-32





SIMPLE ELBOW
 BASİT DİRSEKLİK
 ПРОСТОЙ ЛОКОТЬ

OST 1005
OD 203



It is designed as suitable for use on both sides, by combining pressure and heat that helping to pass the feeling of ache and also accelerates the healing process.

INDICATIONS

- Provides protection and support for the elbow
- Providing heat and compression
- Protection against trauma

UNIVERSAL SIZE

S	M	L	XL	2XL
22-24	24-26	26-28	28-30	30-32



SPORTY ELASTIC KNITTED ELBOW
 SPORTİF ELASTİK ÖRME DİRSEKLİK **OD 204**




It is manufactured from respirable and elastic woven fabric. It provides local compression so reduce tension and pain on tendon.

INDICATIONS

- Intra articular liquid increase
- Provides protection and support for the elbow
- Arthritis
- Support gives pain and stress relief for everyday use

S	M	L	XL
22-25	25-28	28-31	31-34





KNITTED ELBOW

ÖRME DİRSEKLİK

ТРИКОТАЖНАЯ
ПОДСТАВКА ПОД ЛОКОТЬ

OD 205



Manufactured from respirable and elastic woven fabric and designed anatomically. Local compression reduces tendon tension and so the pain.

INDICATIONS

- Intra articular liquid increase
- Provides protection and support for the elbow
- Arthritis
- Support gives pain and stress relief for everyday use

S	M	L	XL	2XL
22-24	24-26	26-28	28-30	30-32



PADDED SPORTS ELBOW BRACE

PEDLİ SPORCU DİRSEKLİK

СПОРТИВНЫЙ ЛОКТЕВОЙ
БАНДАЖ С ПОДУШКОЙ

OD 206



It is designed as suitable for use on both sides. With special pad on elbow protects against sudden impacts. It provides compression and also temperature on elbow to.

INDICATIONS

- Protection during and after the sportive activities
- Protection against trauma

S	M	L	XL	2XL
22-24	24-26	26-28	28-30	30-32





ELBOW IMMOBILISER

DİRSEK İMMOBİLİZER
ИММОБИЛИЗАТОР ЛОКТЯ

OD 207



It provides immobilization extension of elbow. It has removable aluminum supports that provide stabilization in the posterior, medial and lateral. It is supportive for elbow injuries and after surgery. It is universal size and available in 20 cm, 25 cm and 30 cm lengths.

INDICATIONS

- Post-trauma, recovery and preservation purposes
- Elbow joint injuries

20 cm

20 cm

20 cm



CUBITAL TUNNEL ELBOW SUPPORT

KUBİTAL TÜNEL DİRSEKLİK
КУБИТАЛЬНЫЙ
ТУННЕЛЬ КОЛЕНА ОПОРА

OST 1007



There is a 20 degree inclined steel support inside in the elbow. It prevents the elbow joint and ulnar nerve from bending, thus relieves the symptoms. If the portable pad inner elbow is used, protect for impacts.

INDICATIONS

- In the cubital nerve compression
- Rehabilitation after elbow surgery

STD 1 (S/M)

22-26

STD 2 (L/XL/2XL)

26-32





ARM STRAP PRACTICAL

KOL ASKISI PRATİK

РЕМЕНЬ НА РУКУ
ПРАКТИЧЕСКИЙ

OST 1009



With practical designed provides patient's arm to be fixed parallel to the near trunk. It can be used for both sides.

INDICATIONS

- Post-operative rehabilitation
- In rehabilitation of forearm diseases
- Post-trauma, recovery and preservation purposes



UNIVERSAL SIZE



ARM STRAP MESH

KOL ASKISI FİLELİ

НАРУЧНЫЙ
РЕМЕНЬ С СЕТКОЙ

OD 208



It is manufactured from perforated fabric. Product's design provides patient's arm to be fixed parallel to the near trunk. It can be used for both sides.

INDICATIONS

- Post-operative rehabilitation
- In rehabilitation of forearm diseases
- Post-trauma, recovery and preservation purposes



S	M	L	XL	2XL
30-33	33-36	36-39	39-42	42-45



ARM STRAP WITH MESH WAIST BELT

KOL ASKISI FİLELİ BELDEN KEMERLİ
 РЕМЕНЬ НА РУКУ С СЕТЧАТЫМ ПОЯСOM

OD 209



It is manufactured from perforated fabric. Product's design provides patient's arm to be fixed parallel to the near trunk. It can use for both sides. There is a belt for daily use.

INDICATIONS

- Post-operative rehabilitation
- In rehabilitation of forearm diseases
- Post-trauma, recovery and preservation purposes

S	M	L	XL	2XL
30-33	33-36	36-39	39-42	42-45



ARM STRAP LUXURY

KOL ASKISI LÜKS
 РОСКОШНЫЙ НАРУКАВНЫЙ РЕМЕНЬ

OD 210



It is manufactured from sponge laminated fabric. Product's design provides patient's arm to be fixed parallel to the near trunk. It can use for both sides. There is a belt for daily use.

INDICATIONS

- Post-operative rehabilitation
- In rehabilitation of forearm diseases
- Post-trauma, recovery and preservation purposes

S	M	L	XL	2XL
30-33	33-36	36-39	39-42	42-45





30 DEGREE ARM HANGER

30 DERECE KOL ASKISI
ПОДВЕСКА НА 30 ГРАДУСОВ

OST 1010



It is arm sling at 30 degrees and manufactured from sponge laminated fabric. Product 's design provides patient's arm to be fixed parallel to the near trunk. It can use for both sides. There is a belt for daily use.

INDICATIONS

- Rotator cuff in muscle tears
- In shoulder calcification
- Rehabilitation after shoulder implant prosthesis application
- Rehabilitation after shoulder arthroscopy

UNIVERSAL SIZE



30 DEGREE ARM HANGER

30 DERECE KOL ASKISI
ПОДВЕСКА НА 30 ГРАДУСОВ

OST 1011



t is arm sling at 30 degrees and manufactured from sponge laminated fabric. Product 's design provides patient's arm to be fixed parallel to the near trunk. It can use for both sides. There is a belt for daily use.

INDICATIONS

- Rotator cuff in muscle tears
- In shoulder calcification
- Rehabilitation after shoulder implant prosthesis application
- Rehabilitation after shoulder arthroscopy

UNIVERSAL SIZE





45 / 60 DEGREE ARM HANGER

45 / 60 DERECE KOL ASKISI
 ПОДВЕСКА 45/60 ГРАДУСОВ

OST 1012



It is arm sling at 45/60 degrees and manufactured from sponge laminated fabric. Product 's design provides patient's arm to be fixed parallel to the near trunk. It can use for both sides. There is a belt for daily use

INDICATIONS

- Rotator cuff in muscle tears
- In shoulder calcification
- Rehabilitation after shoulder implant prosthesis application
- Rehabilitation after shoulder arthroscopy

UNIVERSAL SIZE



STIMSON BANDAGE

STIMSON BANDAJI

АКРОМИАЛЬНО-КЛЮЧИЧНЫЙ
 БАНДАЖ (БАНДАЖ СТИМСОНА)

OD 211



It provides stabilization between acromio and clavicle. It is generally used to stabilize after humerus operations. It is manufactured from terry cloth. The shoulder pad provides a balanced distribution of the load. Rotating buckle system allows comfortable use of the product.

INDICATIONS

- Acromio-Clavicle dislocations and subluxations
- Post-trauma, recovery and preservation purposes

S **M** **L**
 30-35 35-40 40-45





VELPO BANDAGE

VELPO BANDAJI
БАНДАЖ ВЕЛПО

OD 213
OST 1019



It is made of sponge fabric, provides immobilization on shoulder and upper arm adduction. The bandage promotes stabilization by encircling around the chest. Velcro on the front provides comfortable usage.

INDICATIONS

- Shoulder dislocation – subluxation
- In the tears of rotator cuff muscles
- In humerus (arm bone) and scapula (scapula) fractures
- Injury to the soft tissue of the shoulder and arm

UNIVERSAL SIZE

STD 1 (S/M) **30-36**

STD 2 (L/XL/2XL) **36-45**

S	M	L	XL	2XL
30-33	33-36	36-39	39-42	42-45



SHOULDER DETECTION BANDAGE

OMUZ TESPİT BANDAJI
БАНДАЖ ДЛЯ ОБНАРУЖЕНИЯ ПЛЕЧ

OST 1014



It is a shoulder bandage that protects the heat and improves blood circulation through the compression applied in the shoulder area. It prevents edema due to the applied pressure.

INDICATIONS

- Shoulder bones fractures
- In shoulder calcification
- Rehabilitation after shoulder implant prosthesis application
- In need of hot compress to shoulder

UNIVERSAL SIZE





CLAVICLE BANDAGE

KLAVÍKULA BANDAŽI
БАНДАЖ НА КЛЮЧИЦУ

OST 1017



It is made of sponge fabric with a terry inner side. It's enable for shoulders to be pulled back and smooth posture with velcro adjustment. Clavicle bandage gives accurate immobilization and prevents dislocation by its construction.

INDICATIONS

- Clavicle fractures
- Clavicle and sterno-clavicle dislocations and subluxations
- Used in shoulder joint injuries and posture disorders (postural kyphosis)



STD 1 (S/M) 26-41

STD 2 (L/XL/2XL) 41-61



CLAVICLE BANDAGE WITH PLASTIC SUPPORT

KLAVÍKULA BANDAŽI
PLASTÍK DESTEKLİ

OST 1018

БАНДАЖ НА КЛЮЧИЦУ С ПЛАСТИКОВОЙ ПОДДЕРЖКОЙ



It is made of sponge fabric with a terry inner side. Include extra plastic support. It's enable for shoulders to be pulled back and smooth posture with velcro adjustment. Clavicle bandage gives accurate immobilization and prevents dislocation by its construction.

INDICATIONS

- Clavicle fractures
- Clavicle and sterno-clavicle dislocations and subluxations
- Used in shoulder joint injuries and posture disorders (postural kyphosis)



STD 1 (S/M) 26-41

STD 2 (L/XL/2XL) 41-61



ARM ABDUCTION ORTHOSIS PLUS

KOL ABDÜKSİYON ORTEZİ PLUS
ОРТОЗ ОТВЕДЕНИЯ РУКИ ПЛЮС

OST 1013



To stabilize the shoulder after surgery, a single device immobilizes the shoulder in various abduction positions.

INDICATIONS

- Immobilizes the shoulder at 0°, 15°, 30°, 60°, or 90° abduction.
- Optional External Rotation
Accessory locks the shoulder at 10°, 25°, and 30° external rotation.
- Open design increases airflow and patient comfort.

UNIVERSAL SIZE



HUMERUS BRACE

HUMERUS BREYSİ
ПЛЕЧЕВАЯ КОСТЬ

OD 212



It is made of sponge fabric with a terry inner side. It is a shoulder orthosis covered with polythene material to keep humerus immobilized. Used for dislocations, protection and stabilization of the humerus after surgery and plaster encasements. Has easy usage and application.

INDICATIONS

- Protection and stabilization in shoulder dislocations
- In humeral head fractures
- Rehabilitation after surgical intervention
- Rehabilitation after shoulder plaster application

S M L
22 - 26 26 - 30 30 - 34

LEFT - RIGHT





ELBOW CONTRACTOR ORTHOSIS

DİRSEK KONTRAKTÖR ORTEZİ

ЛОКТОВЫЙ КОНТРАКТОР ОРТОЗ (ДЕРЖАТЬ)

OST 1020



It is a degree adjustable elbow orthosis. Through the mechanism in the elbow area, the desired adjustment can be made. Flexion (0-15-30-45-60-75-90-105) and extension (-15-0-15-30-45-60-75-90) can be adjusted in degrees. Depending on the orthosis, length can be adjusted. Shoulder strap provides comfort in use

INDICATIONS

- In internal or external epicondyle fractures
- After ligament injuries
- Rehabilitation after elbow implant prosthesis application
- Conservative treatment of elbow dislocation

STD 1 (S/M) **30-36**

STD 2 (L/XL/2XL) **36-45**



HUMERUS BRACE PLUS

HUMERUS BREYSİ PLUS

OST 1022



It is a shoulder orthosis made of sponge fabric with a terry inner side. The plastic and aluminum supports inside it keep the humerus immobilized. These supports can be applied to the desired shoulder size. Used for dislocations, protection and stabilization of the humerus after surgery and plaster encasements. It has an optional arm sling.

INDICATIONS

- Protection and stabilization in shoulder dislocations
- In humeral head fractures
- Rehabilitation after surgical intervention
- Rehabilitation after shoulder plaster application

UNIVERSAL SIZE

