



ORFACE[®]
ORTHOPEDIC PRODUCTS

**FINGER SPLINTS
AND WRIST SPLINTS**



ALUMINUM STRAIGHT SPLINT

ALÜMİNYUM DÜZ ATEL
 АЛЮМИНИЕВАЯ ПРЯМАЯ ШИНА

OST 1023




Manufactured from special aluminium with inner plastazote coating. Provides finger immobilization in extension

INDICATIONS

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma

UNIVERSAL SIZE



BASEBALL SPLINT

BEYZBOL ATEL
 БЕЙСБОЛ АТХЕЛ

OPE 100



It is a finger splint and manufactured from special aluminium with inner plastazote coating. Provides finger immobilization in extension.

INDICATIONS

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma

S
4,5-5

M
5-6

L
6-7





BASEBALL FINGER SPLINT

BEZBOL PARMAK ATELİ
БЕЙБОЛЬНЫЙ РАЗДЕЛ ПАЛЬЦЕВ

OST 1024



Made of polyester fabric and under there is an aluminium atele under the finger to keep immobilize. It can be adjustable for finger with velcro on it. It can use for all fingers

INDICATIONS

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma

UNIVERSAL SIZE



FROG SPLINT

KURBAĞA ATEL
ФРОГ СПЛИНТ

OPE 101



It is a finger splint and manufactured from special aluminium with inner plastazote coating. Provides finger immobilization in extension.

INDICATIONS

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma

S	M	L
4,5-5	5-6	6-7





STABILISATION SPLINT

STABİLİZASYON ATEL

СТАБИЛИЗАЦИОННАЯ ШИНА

OPE 102




It is a finger splint and manufactured from special aluminium with inner plastazote coating. Provides finger immobilization in extension.

INDICATIONS

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma

S	M	L
4,5-5	5-6	6-7



ROTECTION FINGER SPLINT

KORUMA PARMAK ATELİ

ЗАЩИТНАЯ ШИНА ДЛЯ ПАЛЬЦЕВ

OPE 123




It is a finger splint and manufactured from special aluminium with inner plastazote coating. Provides finger immobilization in extension.

INDICATIONS

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma

S	M	L
4,5-5	5-6	6-7





CURVED FINGER SPLINT

KAVISLI PARMAK ATELİ
 ИЗОГНУТАЯ ШИНА ДЛЯ ПАЛЬЦЕВ

OPE 124




is a finger splint and manufactured from special aluminium with inner plastazote coating. Provides finger immobilization in extension.

INDICATIONS

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma

S	M	L
4,5-5	5-6	6-7



MALLET FINGER

MALLET FINGER
 ПАЛЬЧИК МАЛЛЕТА

OPE 103



It is a finger orthosis manufactured from special plastic material. The finger is stable in the DIP joint extension, mobile from PIP joint on.

INDICATIONS

- Mallet finger deformity
- PIP joint fractures, traumas and sprains

1	2	3	4
4,5-5	4,5-5	5-5,5	5,5-6
5	5,5	6	7
6-6,5	7,5-8	8,8,5	8,5-9





3 POINT FINGER SPLINT

3 NOKTA PARMAK ATELİ
3-X ТОЧЕЧНАЯ ШИНА ДЛЯ ПАЛЬЦЕВ

OPE 104



It is made of plastic, it's help to stabilizes and aligns the PIP and DIP joints.

INDICATIONS

- Hammer finger
- Butonary disorders
- Fractures
- Arthritis

1-2-3	3-4-5	5-6-7
4,3-4,9	4,9-5,7	5,7-6,2
7-8-9	9-10-11	12-13-14
6,2-6,8	6,8-7,4	7,4-8,2



DOUBLE FINGER BANDAGE

İKİLİ PARMAK BANDAJI
ДВОЙНОЙ БАНДАЖ НА ПАЛЬЦЫ

OST 125



It is an auxiliary bandage consisting of 2 parts that stabilizes 2 finger.

INDICATIONS

- Finger trauma
- Resumption of activity and sport

UNIVERSAL SIZE






FINGER WEIGHT BANDAGE

PARMAK AĞIRLIK BANDAJI
 БАНДАЖ НА ПАЛЬЦЫ

OST 1026



UNIVERSAL SIZE



BUTONIEL SPLINT

BUTONİYER ATEL
 ШИНА ИЗ БУТОНИРА

OPE 106




Finger splint manufactured from special plastic material, keeping the finger in extension according to the 3 point principle.

INDICATIONS

- To correct finger muscle immobility
- Rehabilitation after surgical intervention

S	M	L
4,5-5	5-6	6-7





SPORT FINGER BANDAGE
SPORTİF PARMAK BANDAĞI
 СПОРТИВНЫЙ БАНДАЖ ДЛЯ ПАЛЬЦЕВ

OPE 105




It is an auxiliary bandage consisting of 2 parts that stabilizes 2 finger.

INDICATIONS

- During activity and sport

1 4,5-5	2 4,5-5	3 5,5-5	4 5,5-6	5 6-6,5
6 7,5-8	7 8-8,5	8 8,5-9	9 9-9,5	10 9,5-10



OPE 107

DYNAMIC FINGER REHABILITATION SPLINT
DİNAMİK PARMAK REHABİLİTASYON ATELİ
 ДИНАМИЧЕСКАЯ ШИНА ДЛЯ РЕАБИЛИТАЦИИ ПАЛЬЦЕВ



It is a dynamic splint like coil design for PIP and DIP joint flexion contractures deformity. Active flexion maintains assisted extension in the PIP joint while allowing and affecting the PIP joint with a possible flexion deformity can be used for soft tissue injury.

INDICATIONS

- For PIP and DIP joint flexion
- In injuries related to PIP and DIP
- Rehabilitation after surgical intervention

S 4,5-5	M 5-6	L 6-7
------------	----------	----------





OPE 108

DOUBLE DYNAMIC FINGER REHABILITATION SPLINT

İKİLİ DİNAMİK PARMAK REHABİLİTASYON ATELİ

ДВОЙНАЯ ДИНАМИЧЕСКАЯ
ШИНА ДЛЯ РЕАБИЛИТАЦИИ ПАЛЬЦЕВ



It is a dynamic splint like coil design for PIP and DIP joint flexion contractures deformity. Active flexion maintains assisted extension in the PIP joint while allowing and affecting the PIP joint with a possible flexion deformity can be used for soft tissue injury.

INDICATIONS

- For PIP and DIP joint flexion
- In injuries related to PIP and DIP
- Rehabilitation after surgical intervention

S	M	L
4,5-5	5-6	6-7



OPE 109

TRIPLE DYNAMIC FINGER REHABILITATION SPLINT

ÜÇLÜ DİNAMİK PARMAK REHABİLİTASYON ATELİ

ТРОЙНАЯ ДИНАМИЧЕСКАЯ
ШИНА ДЛЯ РЕАБИЛИТАЦИИ ПАЛЬЦЕВ



It is a dynamic splint like coil design for PIP and DIP joint flexion contractures deformity. Active flexion maintains assisted extension in the PIP joint while allowing and affecting the PIP joint with a possible flexion deformity can be used for soft tissue injury.

INDICATIONS

- For PIP and DIP joint flexion
- In injuries related to PIP and DIP
- Rehabilitation after surgical intervention

S	M	L
4,5-5	5-6	6-7





OST 1028

MIDDLE FINGER FRACTURE SPLINT

ORTA PARMAK KIRIĞI ATELİ
 ПЕРЕЛОМ СРЕДНЕГО ПАЛЬЦА



Allowing stabilization on wrist and middle finger. It can adjust compression by velcro. There is aluminium atele for maximum stabilization and immobilization inner side

INDICATIONS

- Finger fractures
- Tendon and nerve injuries
- Finger injuries, sprains and trauma
- Rehabilitation after surgical intervention

UNIVERSAL SIZE



TRIGGER FINGER SPLINT

TETİK PARMAK ATELİ
 Шина для триггерного пальца

OST 1029



It's provide stabilization of the finger. Can be adjusted to your desired tightness with the high quality velcro straps on it. There is an aluminum splint inside for maximum stabilization and immobilization.

INDICATIONS

- Finger fractures
- Tendon and nerve damage
- Sprains and injuries to fingers
- Rehabilitation after surgery

UNIVERSAL SIZE





BOXER FRACTURE SPLINT

BOKSÖR KIRIĞI ATELİ
 ШИНА ДЛЯ ПЕРЕЛОМОВ БОКСЕРА

OPE 110
OST 1027



Allowing stabilization on little and ring fingers. It can adjust compression by velcro. There is aluminium atele for maximum stabilization and immobilization inner side.

Indications



INDICATIONS

- MP joint fractures
- Rehabilitation after surgical intervention

UNIVERSAL SIZE			RIGHT / LEFT	
S	M	L	XL	2XL
14-16	16-18	18-20	20-22	22-24



SIMPLE BRACELET

BASİT BİLEKLİK
 ПРОСТОЙ БРАСЛЕТ

OST 1037



It is a simple wrist bandage. It provides protection and heat to the wrist in case of wrist injuries, sprains and trauma. It improves blood circulation by maintaining the body temperature and thus speeds up the healing process.

INDICATIONS

- In wrist strains
- To prevent sports injuries
- In carpal dislocations

UNIVERSAL SIZE





THUMB CLOCK BRACELET **OST 1039**

BAŞPARMAK GEÇMELİ BİLEKLİK
БРАСЛЕТ-ЧАСЫ



It is a simple wrist bandage that put thumb to inside. It provides protection and heat to the wrist in case of wrist injuries, sprains and trauma. It improves blood circulation by maintaining the body temperature and thus speeds up the healing process.

INDICATIONS

- In wrist strains
- To prevent sports injuries
- In carpal dislocations

UNIVERSAL SIZE



KNITTED HAND WRIST SPLINT **OPE 111**

ÖRME EL BİLEK ATELİ
ТРИКОТАЖНАЯ ШИНА ДЛЯ ЗАПЯСТЬЯ



Manufactured from respirable and elastic woven fabric. A type of splint that provides immobilization in neutral position. Include a plastic support under the wrist.

INDICATIONS

- Carpal Tunnel Syndrome
- In tendon and soft tissue injuries
- Rehabilitation after surgical intervention

S	M	L	XL	2XL
14-16	16-18	18-20	20-22	22-24





HAND WRIST SPLINT
EL BİLEK ATELİ
 ШИНА ДЛЯ ЗАПЯСТЬЯ РУКИ

OPE 112



Features an anatomically shaped aluminum support that immobilizes the wrist in a neutral position. Includes an additional bandage that enhances stabilization around the wrist.

INDICATIONS

- Carpal Tunnel Syndrome
- In tendon and soft tissue injuries
- Rehabilitation after surgical intervention.

S	M	L	XL	2XL
14-16	16-18	18-20	20-22	22-24



HAND WRIST SPLINT LONG
EL BİLEK ATELİ UZUN
 РУКА ЗАПЯСТЬЕ ШИНА ДЛИННЫЕ

OPE 114



Features an anatomically shaped aluminum support that provides immobilization of the wrist in a neutral position. Includes an additional bandage that enhances stabilization around the wrist.

INDICATIONS

- Carpal Tunnel Syndrome
- In tendon and soft tissue injuries
- Rehabilitation after surgical intervention.

S	M	L	XL	2XL
14-16	16-18	18-20	20-22	22-24





OST 1044

WRIST SPLIT WITH THUMB SUPPORT

BAŞPARMAK DESTEKLİ EL BİLEK ATELİ

РАЗДЕЛЕНИЕ ЗАПЯСТЬЯ С
ПОДДЕРЖКОЙ ДЛЯ БОЛЬШОГО ПАЛЬЦА



Features an aluminum support under the wrist and around the thumb that provides effective stabilization. Includes an additional bandage that enhances stabilization around the wrist.

INDICATIONS

- Combined cases of carpal tunnel syndrome and De Quervain's tenosynovitis
- De Quervain's tenosynovitis
- Degeneration of the 1st CMC and MCP joints
- Painful conditions of the wrist and thumb after injury or surgery
- Post-fracture rehabilitation

S	M	L	XL	2XL
14-16	16-18	18-20	20-22	22-24



HAND WRIST SPLINT WITH PRACTICAL SYSTEM

EL BİLEK ATELİ PRATİK SİSTEMLİ

ШИНА ДЛЯ ЗАПЯСТЬЯ РУКИ С ПРАКТИЧНОЙ СИСТЕМОЙ

OPE 113



Made of leather and breathable fabric. Features an anatomically shaped aluminum support that immobilizes the wrist in a neutral position. Easily adjustable thanks to its practical fastening system.

INDICATIONS

- Protection in sports injuries
- Compression and support after sprains and injuries
- Pre- and post-operative ligament injuries

LEFT / RIGHT

S	M	L	XL
14-16	16-19	19-22	22-25



**WRIST SPLIT**

EL BİLEK ATELİ

ЗАПЯСТЬЕ РАЗДЕЛЕНИЕ

OST 1030

Features an anatomically shaped aluminum support that immobilizes the wrist in a neutral position. Includes an additional bandage that enhances stabilization around the wrist. Designed in a standard size and can be individually adjusted with velcro straps.

INDICATIONS

- Carpal tunnel syndrome
- Tendon and soft tissue injuries
- Post-operative rehabilitation

UNIVERSAL SIZE**OST 1045****WRIST SPLIT WITH THUMB SUPPORT**

BAŞPARMAK DESTEKLİ EL BİLEK ATELİ

РАЗДЕЛЕНИЕ ЗАПЯСТЬЯ С
ПОДДЕРЖКОЙ ДЛЯ БОЛЬШОГО ПАЛЬЦА

Features an aluminum support under the wrist and around the thumb for effective stabilization. Suitable for use on both hands by repositioning the aluminum support. Includes a velcro strap that enhances stabilization around the wrist.

INDICATIONS

- De Quervain's syndrome
- Carpal tunnel syndrome
- Post-operative rehabilitation
- Wrist and thumb osteoarthritis
- Tendon and soft tissue injuries

LEFT / RIGHT**UNIVERSAL SIZE**



OST 1031

LACED WRIST SPLIT WITH MESH

BAĞCIKLI EL BİLEK ATELİ FİLELİ

РАЗДЕЛЕНИЕ НА ЗАПЯСТЬЕ СО ШНУРОВКОЙ И СЕТКОЙ



Provides immobilization of the wrist in a neutral position. Features an internal aluminum support and two plastic supports on the upper side for enhanced stability. Includes a laced fastening system that ensures effective stabilization around the wrist. Suitable for use on both hands.

INDICATIONS

- Carpal tunnel syndrome
- Tendon and soft tissue injuries
- Post-operative rehabilitation

UNIVERSAL SIZE



OST 1047

LACED WRIST SPLIT WITH THUMB SUPPORT

BAĞCIKLI BAŞPARMAK DESTEKLİ EL BİLEK ATELİ

РАЗДЕЛИТЕЛЬ НА ЗАПЯСТЬЕ СО ШНУРОВКОЙ И ПОДДЕРЖКОЙ
ДЛЯ БОЛЬШОГО ПАЛЬЦА



Provides immobilization of the wrist and thumb in a neutral position. Features an aluminum support under the wrist and alongside the thumb, as well as two plastic supports on the wrist for enhanced stability. Adjustable with a lace-up fastening system for a secure fit.

INDICATIONS

- De Quervain's syndrome
- Carpal tunnel syndrome
- Post-operative rehabilitation
- Wrist and thumb osteoarthritis
- Tendon and soft tissue injuries

LEFT / RIGHT

UNIVERSAL SIZE





LACED WRIST SPLINT

BAĞCIKLI EL BİLEK ATELİ

РАЗДЕЛЕНИЕ НА ЗАПЯСТЬЕ СО ШРУРОВОЙ

OST 1032



Made of perforated fabric. Provides immobilization of the wrist in a neutral position. Features an internal aluminum support and two plastic supports on the upper side for enhanced stability. Includes a lace-up fastening system that ensures effective stabilization around the wrist. Suitable for use on both hands

INDICATIONS

- Carpal tunnel syndrome
- Tendon and soft tissue injuries
- Post-operative rehabilitation

UNIVERSAL SIZE



LACED THUMB SUPPORT WRIST SPLIT WITH MESH

BAĞCIKLI BAŞPARMAK DESTEKLİ
EL BİLEK ATELİ FİLELİ

ПОДДЕРЖКА ДЛЯ БОЛЬШОГО БОЛЬШОГО
пальца со шнуровой. Шина на запястье с сеткой.

OST 1046



Provides immobilization of the wrist and thumb in a neutral position. Features an aluminum support under the wrist and alongside the thumb, along with two plastic supports on the wrist for enhanced stability. Adjustable with a lace-up fastening system for a secure fit.

INDICATIONS

- De Quervain's syndrome
- Carpal tunnel syndrome
- Post-operative rehabilitation
- Wrist and thumb osteoarthritis
- Tendon and soft tissue injuries

LEFT / RIGHT

UNIVERSAL SIZE





OST 1033

WRIST SPLIT PLUS (DOUBLE UNDERWIND)

EL BİLEK ATELİ PLUS (ÇİFT BALENLİ)

WRIST SPLIT PLUS (ДВОЙНОЙ ПОДВЕТР)



Provides immobilization of the wrist in a neutral position. Features aluminum supports on both the upper and lower sides of the wrist for enhanced stability. The aluminum splint can be repositioned, allowing use on both hands

INDICATIONS

- Carpal tunnel syndrome
- Tendon and soft tissue injuries
- Post-operative rehabilitation

STD 1 (S/M) 14-18

STD 2 (L/XL/2XL) 18-24



OST 1048

WRIST SPLIT WITH THUMB SUPPORT PLUS

BAŞPARMAK DESTEKLİ EL BİLEK ATELİ PLUS

РАЗДЕЛЕНИЕ ЗАПЯСТЬЯ С ПОДДЕРЖКОЙ ДЛЯ БОЛЬШОГО ПАЛЬЦА ПЛЮС



Provides immobilization of the wrist and thumb in a neutral position. Features an aluminum support under the wrist and alongside the thumb, along with two plastic supports on the wrist for enhanced stability. Adjustable with velcro straps for a secure fit.

INDICATIONS

- De Quervain's syndrome
- Carpal tunnel syndrome
- Post-operative rehabilitation
- Wrist and thumb osteoarthritis
- Tendon and soft tissue injuries

STD 1 (S/M) 14-18

STD 2 (L/XL/2XL) 18-24





LONG WRIST SPLINT PLUS
 UZUN EL BİLEK ATELİ PLUS
 ДЛИННАЯ шина на запястье ПЛЮС

OST 1034



Provides immobilization of the wrist in a neutral position. Features an aluminum support on the underside of the wrist and a formable support at the fingertips for enhanced stabilization. Suitable for use on both hands.

INDICATIONS

- Carpal tunnel syndrome
- Tendon and soft tissue injuries
- Post-operative rehabilitation

UNIVERSAL SIZE



FOREARM SPLINT
 ÖN KOL ATELİ
 СПЛИНД ПРЕДПЛЕЧЬЯ

OST 1035



Provides immobilization of the wrist in a neutral position with its extended splint design. Designed to extend from the forearm and wrap around the wrist for enhanced support. Features an aluminum support on the underside for improved stability. The aluminum support can be repositioned, allowing use on both hands.

INDICATIONS

- Carpal tunnel syndrome
- Tendon and soft tissue injuries
- Post-operative rehabilitation

UNIVERSAL SIZE





FOREARM SPLINT

ÖN KOL SPLİNTİ
СПЛИНД ПРЕДПЛЕЧЬЯ

OST 1036



Designed as an arm splint for use in forearm fractures and injuries. Includes an aluminum support that provides effective stabilization.

INDICATIONS

- Forearm and wrist traumas
- Post-operative rehabilitation

UNIVERSAL SIZE



OST 1049

FOREARM SPLIT WITH THUMB SUPPORT

BAŞPARMAK DESTEKLİ ÖN KOL ATELİ
РАЗДЕЛЕНИЕ ПРЕДПЛЕЧЬЯ С
ПОДДЕРЖКОЙ ДЛЯ БОЛЬШОГО ПАЛЬЦА



Provides immobilization of the wrist and thumb in a neutral position with its extended splint design. Designed to extend from the forearm and securely support the wrist. Includes an internal aluminum support for enhanced stability. The aluminum support can be repositioned, allowing use on both hands.

INDICATIONS

- De Quervain's syndrome
- Carpal tunnel syndrome
- Post-operative rehabilitation
- Wrist and thumb osteoarthritis
- Tendon and soft tissue injuries

UNIVERSAL SIZE





DE QUERVAIN SPLINT

DE QUERVAİN ATELİ
 ШИНА ДЕ КЕРВЕНА

OPE 116



Provides immobilization of the thumb. Features an internal aluminum splint that ensures effective stabilization of the thumb. Adjustable with velcro straps for a secure fit around the wrist.

INDICATIONS

- Lateral injuries of the thumb joint
- Sprains
- 1st MCP joint irritation (e.g. rhizarthrosis)

LEFT / RIGHT

S	M	L	XL	2XL
14-16	16-18	18-20	20-22	22-24



THUMB DETECTION SPLINT

BAŞPARMAK TESPİT ATELİ
 ОБНАРУЖЕНИЕ БОЛЬШОГО ПАЛЬЦА SPLEND

OST 1040



Provides immobilization of the thumb. Features an internal aluminum splint that ensures effective stabilization of the thumb. Adjustable with velcro straps for a secure fit around the wrist.

INDICATIONS

- CMC and MCP joint irritation
- De Quervain's syndrome
- Ligament injuries

UNIVERSAL SIZE





THUMB DETECTION SPINDLE PLUS

BAŞPARMAK TESPİT ATELİ PLUS

Шпindelъ обнаружения большого пальца плюс

OST 1041



Provides immobilization of the thumb. Features an internal plastic splint that ensures effective stabilization of the thumb. Adjustable with velcro straps for a secure fit around the wrist.

INDICATIONS

- CMC and MCP joint irritation
- De Quervain's syndrome
- Ligament injuries



UNIVERSAL SIZE



THUMB SPLINT

BAŞPARMAK SPLİNTİ

ШИНА БОЛЬШОГО ПАЛЬЦА

OST 1042



Provides immobilization of the thumb. Features an internal aluminum splint that ensures effective stabilization of the thumb. Adjustable with velcro straps for a secure fit around the wrist.

INDICATIONS

- CMC and MCP joint irritation
- De Quervain's syndrome
- Ligament injuries



UNIVERSAL SIZE

**THUMB SPLINT**

BAŞPARMAK TESPİT ATELİ
 ОБНАРУЖЕНИЕ БОЛЬШОГО ПАЛЬЦА

OST 1043

Restricts movement of the thumb and provides effective stabilization. Features a plastic support that ensures thumb stability. Easily adjustable with high-quality velcro straps for a secure fit around the wrist.

INDICATIONS

- CMC and MCP joint irritation
- De Quervain's syndrome
- Ligament injuries

STD 1 (S/M) **14-18**

STD 2 (L/XL/2XL) **18-24**

**DROP WRIST SPLINT**

DÜŞÜK EL BİLEK ATELİ
 РАЗДЕЛЕНИЕ ЗАПЯСТЬЯ

OST 1053

Restricts movement of the thumb and provides effective stabilization. Features a plastic support that ensures thumb stability. Easily adjustable with high-quality velcro straps for a secure fit around the wrist.

INDICATIONS

- CMC and MCP joint irritation
- De Quervain's syndrome
- Ligament injuries

UNIVERSAL SIZE



THERMOPLASTIC HAND RESTING SPLINT

TERMOPLASTİK EL
İSTIRAHAT SPLİNTİ
ТЕРМОПЛАСТИЧНЫЙ УПОР ДЛЯ РУК

OPE 117



INDICATIONS

- Wrist fractures after plaster casting
- Soft tissue injuries requiring rest
- Carpal tunnel syndrome
- Post-operative rehabilitation
- Tendon and soft tissue injuries

LEFT / RIGHT

S	M	L
14-18	18-22	22-26



THERMOPLASTIC HAND WRIST SPLINT

TERMOPLASTİK
EL BİLEK SPLİNTİ
ТЕРМОПЛАСТИЧНАЯ ШИНА ДЛЯ ЗАПЯСТЬЯ РУКИ

OPE 119



Hand and wrist orthosis made of thermoplastic material. Maintains the wrist and fingers in a functional resting position. Provides immobilization of the wrist and is suitable for use during recovery following plaster treatment and in painful wrist conditions.

INDICATIONS

- Wrist fractures after plaster casting
- Soft tissue injuries requiring rest
- Carpal tunnel syndrome
- Post-operative rehabilitation
- Tendon and soft tissue injuries

LEFT / RIGHT

S	M	L
14-18	18-22	22-26





THERMOPLASTIC ANTI-SPASTICITY SPLINT

TERMOPLASTİK
ANTİSPASTİSİTE ATELİ

ТЕРМОПЛАСТИЧНАЯ ШИНА ПРОТИВ СПАСТИКИ

OPE 120



Hand and wrist orthosis made of thermoplastic material with an inner plastazote lining. The grooved finger design and the material's resistance to tension make it suitable for managing hand spasticity.

INDICATIONS

- Spasticity
- Wrist disorders

LEFT / RIGHT

S	M	L
14-18	18-22	22-26



ANTISPASTICITY SPLINT (NEW MODEL)

ANTİSPASTİSİTE ATELİ

АНТИСПАСТИЧЕСКАЯ ШИНА (НОВАЯ МОДЕЛЬ)

OPE 121



Hand and wrist orthosis with a silicone-coated interior. The grooved finger design and tension-resistant structure make it suitable for managing hand spasticity.

INDICATIONS

- Spasticity
- Wrist disorders

LEFT / RIGHT

S	M	L
14-18	18-22	22-26





THERMOPLASTY DYNAMIC HAND SPLINT

TERMOPLASTİ DİNAMİK EL ATELI

ТЕРМОПЛАСТИЧЕСКАЯ
ДИНАМИЧЕСКАЯ ШИНА ДЛЯ РУК

OPE 122



Hand and wrist orthosis made of thermoplastic material. Maintains the wrist and fingers in passive extension while allowing active finger flexion.

INDICATIONS

- Muscle immobility correction
- Ulnar nerve injuries
- Post-operative rehabilitation

LEFT / RIGHT

S	M	L
14-18	18-22	22-26



THERMOPLASTIC KLEINERT SPLINT

TERMOPLASTİK KLEINERT ATELI

ТЕРМОПЛАСТИК КЛЕЙНЕРТ СПЛИ

OST 1052



Hand and wrist orthosis made of thermoplastic material. Allows passive flexion while enabling active extension.

INDICATIONS

- Flexor tendon injuries
- Post-operative rehabilitation

LEFT / RIGHT

S	M
14-18	18-26





FABRIC HAND REST SPLINT

KUMAŞLI EL İSTIRAHAT SPLİNTİ
 ТКАНЕВАЯ ПОДДЕРЖКА ДЛЯ РУК SPLIN

OST 1050



Hand and wrist orthosis that keeps the wrist and fingers in a resting position. Provides stabilization after plaster casting and in painful wrist conditions.

INDICATIONS

- Wrist fractures after plaster casting
- Soft tissue injuries requiring rest
- Carpal tunnel syndrome
- Post-operative rehabilitation
- Tendon and soft tissue injuries

LEFT / RIGHT

UNIVERSAL SIZE



DYNAMIC HAND SPLINT

DİNAMİK EL ATELİ
 ДИНАМИЧЕСКИЙ РУЧНОЙ СПЛИТ

OST 1051



Hand and wrist orthosis made of thermoplastic material. Maintains the wrist and fingers in passive extension while allowing active finger flexion.

INDICATIONS

- Muscle immobility correction
- Ulnar nerve injuries
- Post-operative rehabilitation

UNIVERSAL SIZE

